

DTIC FILE COPY

4

**CHEMICAL
RESEARCH,
DEVELOPMENT &
ENGINEERING
CENTER**

CRDEC-TR-88160

**LIGHT-OFF TEMPERATURE DETERMINATION
OF OXIDATION CATALYST USING
FOURIER TRANSFORM INFRARED TECHNIQUE**

by **Chen C. Hsu, Ph.D.**
RESEARCH DIRECTORATE

August 1988

DTIC
ELECTE
OCT 1 1 1988
S H D

**U.S. ARMY
ARMAMENT
MUNITIONS
CHEMICAL COMMAND**



Aberdeen Proving Ground, Maryland 21010-5423

DEFENSE DOCUMENTATION CENTER

Approved for public release;

distribution is unlimited.

AD-A199 596

88-07-708

Disclaimer

The findings in this report are not to be construed as an official Department of the Army position unless so designated by other authorizing documents.

Distribution Statement

Approved for public release; distribution is unlimited.

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE

REPORT DOCUMENTATION PAGE				Form Approved OMB No. 0704-0188	
1a. REPORT SECURITY CLASSIFICATION UNCLASSIFIED			1b. RESTRICTIVE MARKINGS		
2a. SECURITY CLASSIFICATION AUTHORITY			3. DISTRIBUTION / AVAILABILITY OF REPORT		
2b. DECLASSIFICATION / DOWNGRADING SCHEDULE			Approved for public release; distribution is unlimited.		
4. PERFORMING ORGANIZATION REPORT NUMBER(S) CRDEC-TR-88160			5. MONITORING ORGANIZATION REPORT NUMBER(S)		
6a. NAME OF PERFORMING ORGANIZATION CRDEC		6b. OFFICE SYMBOL (If applicable) SMCCR-RSC-A	7a. NAME OF MONITORING ORGANIZATION		
6c. ADDRESS (City, State, and ZIP Code) Aberdeen Proving Ground, MD 21010-5423			7b. ADDRESS (City, State, and ZIP Code)		
8a. NAME OF FUNDING / SPONSORING ORGANIZATION CRDEC		8b. OFFICE SYMBOL (If applicable) SMCCR-RSC-A	9. PROCUREMENT INSTRUMENT IDENTIFICATION NUMBER		
8c. ADDRESS (City, State, and ZIP Code) Aberdeen Proving Ground, MD 21010-5423			10. SOURCE OF FUNDING NUMBERS		
		PROGRAM ELEMENT NO.	PROJECT NO. 1C161102	TASK NO. A71A	WORK UNIT ACCESSION NO.
11. TITLE (Include Security Classification) Light-Off Temperature Determination of Oxidation Catalyst Using Fourier Transform Infrared Technique					
12. PERSONAL AUTHOR(S) Hsu, Chen C., Ph.D.					
13a. TYPE OF REPORT Technical		13b. TIME COVERED FROM 87 Jul TO 87 Dec		14. DATE OF REPORT (Year, Month, Day) 1988 August	
15. PAGE COUNT 15					
16. SUPPLEMENTARY NOTATION This document was submitted to Journal of Catalysis, May 1988.					
17. COSATI CODES			18. SUBJECT TERMS (Continue on reverse if necessary and identify by block number)		
FIELD	GROUP	SUB-GROUP	Light-off temperature, Oxidation)		
15	06	03	Fourier transform, IR spectroscopy, Catalyst.		
			Oxidation products,		
19. ABSTRACT (Continue on reverse if necessary and identify by block number) Light-off temperatures of oxidation catalysts are considered as one of the important parameters for catalyst performance evaluation. In this study, the in situ Fourier Transform IR technique was developed and used to determine the light-off temperatures and reaction products of three three-way automotive catalysts with 20 torr monomethylamine in air at 0.5 L/min flow rate. Light-off temperatures were 140, 143, and 170 °C for Davison, Allied-Signal, and Degussa oxidation catalysts, respectively. CO, CO ₂ , H ₂ O, and nitric acid were the major oxidation products. The activation energies of formation of CO and CO ₂ on the catalysts were also determined.					
20. DISTRIBUTION / AVAILABILITY OF ABSTRACT <input checked="" type="checkbox"/> UNCLASSIFIED/UNLIMITED <input type="checkbox"/> SAME AS RPT. <input type="checkbox"/> DTIC USERS			21. ABSTRACT SECURITY CLASSIFICATION UNCLASSIFIED		
22a. NAME OF RESPONSIBLE INDIVIDUAL SANDRA J. JOHNSON			22b. TELEPHONE (Include Area Code) (301) 671-2914		22c. OFFICE SYMBOL SMCCR-SPS-T

UNCLASSIFIED



PREFACE

The work described in this report was authorized under Project No. 1C161102A71A, Research in CW/CB Defense. This work was started in July 1987 and completed in December 1987.

The use of trade names or manufacturers' names in this report does not constitute an official endorsement of any commercial products. This report may not be cited for purposes of advertisement.

Reproduction of this document in whole or in part is prohibited except with permission of the Commander, U.S. Army Chemical Research, Development and Engineering Center, ATTN: SMCCR-SPS-T, Aberdeen Proving Ground, Maryland 21010-5423. However, the Defense Technical Information Center and the National Technical Information Service are authorized to reproduce the document for U.S. Government purposes.

This report has been approved for release to the public.

Acknowledgments

The author thanks Marie A. Yeh of C. Milton Wright High School, Bel Air, MD, who performed the experiments during her tenure as a Department of Defense Science and Engineering Apprentice during the summer of 1987. The author especially thanks Dr. G. R. Lester of Allied-Signal, Dr. R. F. Freese of Degussa, and Dr. P. A. Smith of Davison, who provided the three-way automotive catalysts used in this study.

DTIC
COPY
1987

Accession For	
NTIS CRA&I	<input checked="checked" type="checkbox"/>
DTIC TAB	<input type="checkbox"/>
Unannounced	<input type="checkbox"/>
Justification	
By	
Distribution/	
Availability	
Availability	
Dist	
A-1	

Blank

CONTENTS

		Page
1.	INTRODUCTION	7
2.	EXPERIMENTAL PROCEDURES	7
3.	RESULTS AND DISCUSSION	9
4.	SUMMARY	13
	LITERATURE CITED	15

LIST OF FIGURES

Figure

1.	FTIR Catalyst Light-Off Temperature Measurement System	8
2.	Schematic Diagram of the Flow IR Cell	8
3.	IR Spectrum of MMA with Davison Catalyst at 94 °C	10
4.	IR Spectrum of MMA with Davison Catalyst at 190 °C	10
5.	IR Spectrum of MMA with Davison Catalyst at 261 °C	11
6.	Arrhenius Plot of CO on Davison Oxidation Catalyst	12
7.	Arrhenius Plot of CO ₂ on Davison Oxidation Catalyst	12

LIST OF TABLES

Table

1.	Light-Off Temperatures of Three-Way Automotive Catalysts ...	11
2.	Activation Energy of Formation, Kcal/mol	13

Blank

LIGHT-OFF TEMPERATURE DETERMINATION OF OXIDATION CATALYST USING FOURIER TRANSFORM INFRARED TECHNIQUE


1. INTRODUCTION

There are various techniques to evaluate the performance of oxidation catalysts. Measurements of surface composition, surface structure, surface area, porosity, acidity, and dispersion of active metals of oxidation catalysts, etc., can be related to the catalyst performance.¹ Other methods such as the decoloration of the indigo carmine solution to correlate the catalytic activity with the flash point² and the recorder-equipped calorimeter technique to determine the ignition time, maximum temperature, and maximum combustion temperature³ were used. For an oxidation catalyst, however, one direct measurement of its performance is to determine its light-off temperature, the temperature at which significant oxidation reactions occur. In general, it is true that the lower the light-off temperature, the more effective will be the catalyst performance.⁴⁻⁵ Indeed, this correlation was used in the past to select the optimum reduction temperature and the duration of reduction for the preparation of the best Ni-containing catalyst.⁶ We developed an in situ Fourier Transform IR (FTIR) technique for simultaneously determining the light-off temperatures and identifying catalytic oxidation products of oxidation catalysts. Using this newly developed technique, three three-way automotive catalysts were evaluated. The details of the subject technique and the results of catalyst performance evaluation are described below.

2. EXPERIMENTAL PROCEDURES


Anhydrous monomethylamine (MMA) used in this study was supplied by Matheson Gas Products with a purity of 99.99%. MMA was used as received without further purification. Three three-way automotive oxidation catalysts were provided by Allied-Signal (UOP# 4576-136), Davison (alumina catalysts code 701), and Degussa (TWC-1, MS-599). Bead-type catalyst samples were dried overnight in an oven at 110 °C before use.

The schematic diagram of the in situ FTIR catalyst light-off temperature measurement system is shown in Figure 1. The measurement system consists of a flow IR cell, a K-type thermocouple with a chart recorder and digital thermometer display, a mass flow controller, a lecture bottle of MMA, a MKS pressure gauge, and Nicolet FTIR Model 60SX. The details of the flow IR cell are illustrated in Figure 2. This flow IR cell has 10 cm optical pathlength and a 25 x 4 mm NaCl window, the heating element and sample holder, the gas inlet and outlet, and the thermocouple. The thermocouple is positioned in the proximity of the sample to measure the temperature of an oxidation reaction. For accurate determination of the light-off temperature of a catalyst, a chart recorder was used along with a digital thermometer display.



- | | | |
|-------------------------|--------------------|--|
| 1. Compressed Air | 5. Exhaust | 9. Thermocouple |
| 2. Mass Flow Controller | 6. IR Cell | 10. Sample Compartment of
Nicolet FTIR Model 60SX |
| 3. Gas Lecture Bottle | 7. IR Window | 11. Chart Recorder |
| 4. MKS Pressure Gauge | 8. Catalyst Sample | 12. Digital Thermometer
Display |

Figure 1. FTIR Catalyst Light-Off Temperature Measurement System



- | | |
|-----------------|--------------------------------------|
| 1. Gas Inlet | 5. Catalyst Sample |
| 2. Gas Outlet | 6. Electrical Leads for Heating |
| 3. IR Window | 7. Heating Element and Sample Holder |
| 4. Thermocouple | 8. IR Beam |

Figure 2. Schematic Diagram of the Flow IR Cell

Initially, a correlation between the temperature setting on the heater controller and the temperature reading on the digital thermocouple display was obtained. The temperature reading was about 30 to 40% of the temperature setting on the heater controller. To run an experiment, the IR background spectra were taken at various temperatures with the flow IR cell purged with dried air at 0.5 L/min flow rate, which was automatically controlled by a Matheson mass flow controller. Background spectra were taken in the mixture of dried air and 20 torr MMA in the absence of catalyst. The background spectra essentially remain the same in the temperature range of 25-363 °C indicative of no catalytic activity due to the heating element. For catalyst evaluation, three catalyst beads were placed on the sample holder (the heating element with a V-shape), and the thermocouple was positioned in the proximity of catalyst beads. After the light-off temperature measurement system was assembled, the flow IR cell was evacuated with a vacuum pump, and MMA was injected and maintained at 20 torr. The system was disconnected with the vacuum pump, and the IR cell was brought to ambient pressure with compressed air at 0.5 L/min flow rate. The heater controller was set at around 200 °C (equivalent to the thermocouple reading of 80 °C), and the setting was increased 10° every 15-20 min. During this time, an IR spectrum was taken between 4,000 and 400 cm^{-1} with 4 cm^{-1} resolution. Normally, 32 scans were taken for each IR spectrum at every temperature setting. Significant oxidation reaction was detected around 140 °C thermocouple reading on the digital thermometer display. The determination of light-off temperature was carried out using the chart of the temperature recorder. The light-off temperatures were reproducible with repeated on and off cycles of MMA. These procedures were used for three three-way automotive oxidation catalysts provided by Davison, Allied-Signal, and Degussa.

3. RESULTS AND DISCUSSION

Typical IR spectra at various temperatures are illustrated in Figures 3-5 with Davison three-way catalyst samples. At 94 °C, the IR features are very much the same as MMA as shown in Figure 3. When the temperature was increased to 190 °C (as shown in Figure 4), peaks around 2,154 and 2,351 cm^{-1} respectively corresponding to CO and CO₂, and peaks around 1,540 and 3,770 cm^{-1} corresponding to H₂O grow substantially. As the temperature reaches 261 °C, the formation of HNO₃ and the effect of water on nitric acid⁷⁻⁹ corresponding to the bands centered around 2,964, 1,700, 1,320, 760 cm^{-1} appear to be the major reaction products, Figure 5. For each sample, at least three experiments were conducted, and the light-off temperatures were determined using the charts from the temperature recorder. The average light-off temperatures for three three-way automotive catalysts are shown in Table 1. In the increasing light-off temperature, the order is Davison, Allied-Signal, and Degussa catalysts. For Davison and Allied-Signal samples, the light-off temperatures are very close to each other. However, for the Degussa sample, the light-off temperature is about 30 °C higher than that of Davison and Allied-Signal.




Figure 3. IR Spectrum of MMA with Davison Catalyst at 94 °C




Figure 4. IR Spectrum of MMA with Davison Catalyst at 190 °C




Figure 5. IR Spectrum of MMA with Davison Catalyst at 261 °C

Table 1. Light-Off Temperatures of Three-Way Automotive Catalysts

Davison	140 °C
Allied-Signal	143 °C
Degussa	170 °C

The absorbance of CO and CO₂ were taken from the IR spectra at various temperatures, and Arrhenius plots were obtained as shown in Figures 6-7. The activation energies of formation of CO and CO₂ were found using the relation of slope = $-E_a/R$, where E_a is the activation energy and R is the ideal gas constant, which is equal to 1.987 cal/mol/K. Table 2 shows the activation energies of formation of CO and CO₂ for three three-way automotive catalysts. In general, the activation energy of formation of CO₂ is lower than that of CO, which is consistent with the free energy of formation of CO₂ which is lower than that of CO. Note that the lower the light-off temperature, the lower the activation energy of formation for both CO and CO₂ from the three-way automotive catalysts investigated.




Figure 6. Arrhenius Plot of CO on Davison Oxidation Catalyst




Figure 7. Arrhenius Plot of CO₂ on Davison Oxidation Catalyst

Table 2. Activation Energy of Formation, Kcal/mol

	CO ₂	CO
Davison	2.29	6.41
Allied-Signal	2.68	12.07
Degussa	5.76	21.60

4. SUMMARY

The results are summarized as follows:

- An in situ FTIR technique was developed and demonstrated for the light-off temperature determination and the identification of reaction products of oxidation catalysts.

- Light-off temperatures with 20 torr MMA for Davison, Allied-Signal, and Degussa three-way automotive catalysts were 140, 143, and 170 °C in air with 0.5 L/min flow rate, respectively.

- The reaction products were CO, CO₂, H₂O, and HNO₃ for three three-way automotive catalysts.

- The activation energies of formation for CO and CO₂ on three three-way catalysts were determined. For CO, they are 6.41, 12.07, and 21.60 Kcal/mol; and for CO₂, they are 2.29, 2.68, and 5.76 Kcal/mol, for Davison, Allied-Signal, and Degussa three-way automotive catalysts, respectively.

Blank

LITERATURE CITED

1. Delannay, F., and Delmon, B., "Methods of Catalyst Characterization: An Overview," In Characterization of Heterogeneous Catalysts, F. Delannay, ed., Marcel Dekker, Inc., New York and Basel, 1984.
2. Krause, A., and Nowakowski, E., "Catalytic Properties and Flash Points of Some Aromatic Hydrocarbons," Chem.-Zig., Chem. App. Vol. 92(21), pp 774-775 (1968).
3. Guveev, A.A., Dovlator, I.A., Kazaryan, S.A., and Naftulin, I.S., "Determination of the Octane Number of Gasoline from Catalytic Reforming According to the Characteristics of Low-Temperature Gas-Phase Oxidation," Neftepererab. Neftekhim. (Moscow) Vol. 9, pp 5-8 (1987).
4. Pfefferle, L.D., and Pfefferle, W.C., "Catalysis in Combustion," Catal. Rev.-Sci. Eng. Vol. 29(2&3), pp 219-267 (1987).
5. Prasad, R., Kennedy, L.A., and Ruckenstein, E., "Catalytic Combustion," Catal. Rev.-Sci. Eng. Vol. 26(1), pp 1-58 (1984).
6. Anderlohr, A., and Hedden, K., "Characterization of Catalyst Activity by Determination of Ignition Temperature of Exothermic Reactions," Chem.-Ing.-Tech. Vol. 51(8), pp 811 (1979).
7. McGraw, G.E., Bernitt, D.L., and Hisatsune, I.C., "Vibrational Spectra of Isotopic Nitric Acid," J. Chem. Phys. Vol. 42(1), pp 237-244 (1965).
8. Ingold, C.K., and Millen, D.J., J. Chem. Soc. Vol. 2612 (1950).
9. Stern, S.A., Mullhaupt, J.T., and Kay, W.B., Chem. Revs. Vol. 60, p 185 (1960).